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Operating Instructions

Mode d'emploi

- » Handschuhboxe Pure mit Schleuse
- » Glove Box Pure with Transfer Chamber

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SICCO Glove Box Pure, Cat. No.: V 1740-08

OPERATING INSTRUCTIONS

1. SUPPLIED WITH	16
2 . SAFETY ADVICE	16
3. TAKING INTO OPERATION	17
4 FUNCTIONALITY	23
5. SERVICE / MAINTENANCE	25
6. DISTURBANCES / REPAIR	26
7. WASTE DISPOSAL	27
8. SPECIFICATIONS	27

SICCO Glove Box Pure

Please read the information listed below thoroughly. We urgently ask you to observe the instructions about safety, use and maintenance.

SICCO Glove Boxes provide safe handling of toxic materials, allow production of chemical compounds and protect sensitive substances from the environment. The transfer chamber allows to introduce products without noticeable change of the atmosphere inside.

The integrated UV lamp works without creating ozone and emits wavelengths of 254 nm to sterilise the surfaces in the workspace.

All panels are made of high-quality and transparent acrylic glass (PMMA). The aluminium profiles are joint together with connectors made of polyamide.

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- 1 SICCO Glove Box Pure with integrated UV lamp and control system to monitor lifetime
- 1 Pair of gloves size 9,75 made of UV resistant EPDM
- 2 connecting kits for gas-filling for SICCO glove boxes with transfer chamber
- 1 Operating Instructions

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- » Do not take a damaged glove box into operation.
- » Check before every use whether the gloves have been damaged.
- » Do not work with opened swing-up door. Just activate the UV lamp with closed swing-up door.
- » Glove Box shall be taken into operation only in rooms with sufficient ventilation.
- » Place the glove box in a reachable working height.
- » Do not expose the glove box to direct sunlight.
- » Do not store any solvents in the glove box.
- » Do not use any sharp or pointed tools inside the glove box.
- » Do not put any explosive or inflammable materials inside the glove box.
- » Do not put any radioactive products inside the glove box.
- » Glove Boxes are not suitable for overpressure or vacuum.
- » The ambient temperature range of $-20\text{ }^{\circ}\text{C}$ to $+70\text{ }^{\circ}\text{C}$ shall not be exceeded.
- » These operating instructions cannot reflect all possible dangers which might occur during use. Users should check and observe all safety regulations.
- » This device must not be used by persons with limited physical, sensory or mental capacities or without experience and knowledge.
- » The UV lamp is a technical work tool and only intended for the use in our Glove Boxes.
- » Do not look directly into the UV lamp when it is switched on.
- » Protect your eyes and skin from UV light! UV-C radiation creates sunburn and painful conjunctivitis within a short time.
- » The unintended use of the device or damage of the casing can cause leakage of harmful UV-C radiation.
- » Devices with obvious damages must not be taken into operation.
- » Make sure to install all connecting kits for gas-filling.
- » Always use the smaller tubing as your tubing for gas-filling and the larger tubing as your aeration tubing.
- » The pressure for lead-in gas shall not exceed 0,3 bar.
- » Connect the aeration tubing at first, then the tubing for gas-filling. Disconnection has to be made in reverse order.
- » The cross section of the tubing must not be constricted.
- » Neither put the tubing for gas-filling nor the aeration tubing under tensile or compressive stress.
- » Avoid unanticipated reactions which might be caused by gas residues upon opening the swing-up door or the transfer chamber.

Attention: Areas which can only be opened with tools, belong to maintenance areas. In case of unauthorized opening, dangerous risks can arise for the user!

Sicco Glove Box Pure

Taking into operation

Remove the packing and the protection cover and place the glove box on a horizontal, plain bench.

Attention: Never lift the glove box at the transfer chamber. Only lift the glove box at the bottom (see picture below).

Just lift here.

Just lift here.

Never lift up here.


Swing-updoor

» Always pull both handles to open the swing-up door. Swing the door panel up to the highest point until it stops.

Please note: Due to the strong magnetic lock it is harder to open the door the longer it was closed before. This is normal and can be compensated by pulling stronger than in everyday use.

» **Attention:** The door stop defines the maximum opening of the swing-up door. Put the door on the door stop. The adjoining picture shows the maximum possible opening of the door


Swapping with Gas

The substances and components which are stored in the glove box can be swapped with gas. The corresponding bores to lead in the gas are located on the side walls and are factory-made. For technical reasons (packing and transport), the tubing and the connecting pieces are enclosed in delivery but not yet mounted. The assembly is quite easy.

Please note: gas which is lighter than air has to be lead in from above. Gas which is heavier than air from the bottom (density of air = 1,2041 kg/m³). This information is important for further use since the tubing for gas-filling has a different cross section then the aeration tubing.

To better understand, the following pictures shows the different tubing.

» Always use the smaller tubing (3) for gas-filling. Use the bigger tubing (4) for aeration only. This way you can prevent the accruegment of overpressure and therefore any damage to your glove box.

» Open the four cable lead-ins (two cable lead-ins on the back side of the glove box and two on the transfer chamber) by unscrewing the sleeve nuts (1). Remove the sealing disks (2). Please keep the sleeve nuts and the sealing disks safe for re-use.


» Cut all tubing with a knife to the requested length. 5 meters of each tubing are included in delivery (for models with transfer chamber 10 m, respectively 15 m, are included). The picture on the right shows the tubing for gas-filling.


» Connect the tubing for gas-filling by sliding them on the hose connectors. Assemble the aeration tubing accordingly. The picture on the right shows the assembly of the tubing for gas-filling.


» **Tip:** Wet the hose connectors with some water or silicone oil.

» If the gas used for swapping is heavier than air the smaller tubing for gas-filling (3) has to be assembled on the lower cable lead-ins. The lager aeration tubing (4) has to be assembled on the upper cable lead-ins.

» For assembling the tubing please use the included hexagon wrench. Fix the four fastening screws by turning them carefully clockwise. The picture on the right shows the aeration tubing.

» If the gas used for swapping is lighter than air the tubing for gas-filling (3) has to be assembled on the upper cable lead-ins. Corresponding to this the aeration tubing (4) has to be assembled on the lower cable lead-ins.

Attention:


- » Transfer the extracted air in a controlled way, e.g. by using an extractor hood. Please make sure that the airflow of the extractor hood does not effect or ruin the process of gas-filling in our glove box.
- » The cross section of the tubing must not be constricted. Make sure to choose a sufficient bending radius for the tubing when installing them. Do not put any components on the tubing which could damage or crush them.
- » We recommend to purge the transfer chamber and the work space with 0,15 bar or with 0,30 bar. The following table shows you the scavenging time:

Scavenging pressure inside the inlet hose	Flow rate approx.	Scavenging timeinside work space minimum	Scavenging time inside transfer chamber mindestens
0,15 bar	3100 liters per hour	6:00 minutes	1:00 minute
0,30 Bar	6200 liters per hour	3:00 minutes	0:30 0:30 minute

Cable lead-ins to lead in cables or rigid tubing

The cable lead-ins on the backside of the chamber allow to lead in any standard cables including mounted plugs or any rigid tubing. Every cable or rigid tubing with an outer diameter from 0,1 mm up to 10 mm can be assembled regardless of the cross section the cable / tubing has.

» Unscrew the four cylinder screws counter-clockwise. A suitable hexagon wrench is included.


» Remove the supporting ring (5) and the supporting disk with centring collar (6).

» Remove the silicone plug (7) and the sealing foil (8).

» Remove the clampring (9).

» Put your electrical device inside the glove box. Lead your cable through the supporting ring (5). Please pay attention to the supporting ring (5) being properly assembled with the notches in direction to yourself in order to fix the cylinder screws later.


Sicco Glove Box Pure

» Guide the plug from inside the glove box through the housing (10) of the cable lead-in.


» Pull the cable of your electrical device until the required position is reached. Slide the clamping ring (9) on the cable.

» Afterwards slide the silicone plug (7) on the cable as well.

» Then slide the supporting disk with centring collar (6) on the cable with the centring collar in direction to yourself.

» Press the silicone plug (7) into the housing (10). While doing so put pressure to the housing (10) from outside as well.

» Shove the supporting disk with centring collar (6) along the cable. Place the clamp ring with ridge on top of the silicone plug (7). Place the supporting ring (5) on the supporting disk with centring collar (6).

Sicco Glove Box Pure

» Fix the four fastening screws by turning them carefully clockwise.

» The assembly of a rigid tubing is identical to the assembly of a cable. Just follow the abovementioned steps.

Transfer Chamber

The transfer chamber is used to introduce further products or tools during a running process. A contamination of the working space or a loss of the filled-in gas can be prevented.

» Close both doors of the transfer chamber.

» If the inner door of the transfer chamber is open, the interlock prevents accidental opening of the outer door.

» Put the products into the transfer chamber.

» A tray (accessory, Cat. No. V 1861-07) facilitates the removal of the items later.

» Close the outer door.

- » Open the inner door of the transfer chamber and remove the products or tools.

4. Functionality

The working principle of the disinfection with UV light is as follows: Especially short-wave UVC rays have an antiseptic effect. Microorganisms such as bacteria, mould and yeast fungus, primordial animals and viruses are inactivated and destroyed through an appropriately long radiation in a wave range of approx. 264 nm. The exposure time depends on the type of microorganism and takes up to approx. 625 seconds. You have to add a heating time of approx. 30 seconds.

All necessary devices such as balances, syringes, pipettes etc. can be installed in the workspace of the glove box through the large front opening before starting with work. Power supply of these devices is made through the rear side cable-lead ins which allow to lead in cables with O. D. 0,1 to 10 mm as well as the plug.

Please consider when placing your devices that surfaces in shaded areas will not be disinfected.

As soon as the door is closed, the workspace can be sterilised. The control unit is easily accessible on the upper side of the glove box. As UV rays are invisible, a control lamp shows if the UV radiator is turned on. Further control lamps monitor the service life of the illuminant and shows in time when an exchange is necessary. At the end of the service life, the UV rays do no longer attain the requested effect.

The acrylic glass panels completely shield the UV rays so there is no radiation emitted to the environment. During operation, you should always either wear UV-resistant gloves (included in the scope of delivery) or mount the closing lid.

The integrated UV radiator has a lifetime of approx. 10.000 hours and should be switched on and off maximum 5 to 6 times per day. It is an ozone-free source of radiation from which no rays below 220 nm emit.


Dosing groups and Irradiation times

The irradiation times are calculated on the largest distance in the box and maximum of lethal dose of the group.

Microorganism	Irradiation time
Bacillus anthracis / anthrax S. enteritidis, S. typhimurium / salmonellae Corynebacterium diphtheriae / diphtheria Pseudomonas aeruginosa / cocci Eberthella Typhosa / typhoid Neisseria Catarrhalis / cocc Pseudomonas aeruginosa / Verdigris bacterium Shigella dysenteriae / Shigellenruhr Bacillus Tuberculi / tuberculosis Infectus Hepatitis / hepatitis Influenza / flu virus Haemophilus Staphylococcus aureus Bacillus dysenteria Clostridium Tetani / tetanus Measles virus	approx. 133 seconds plus heating time
B. cereus spores Sarcinalutea Polio virus 2 Pneumococci	approx. 624 seconds plus heating time

Meaning of the control unit lamps:

- » Main switch lights up: The UV lamp is switched on.
- » Green LED lights up: The UV lamp is still working, the up-to-now-service life is below 6.000 hours.
- » Yellow LED lights up: The remaining service life is still only 5 % of the maximum running hours. An exchange of the UV lamp illuminant should be prepared.
- » Red LED lights up: The maximum service life has been reached. The UV lamp illuminant has to be exchanged as an efficient sterilisation is no longer ensured.
- » Reset button: After successful exchange of the illuminant, the integrated operating hours counter is reset. A correct service life status is shown through the LED lamps for the new illuminant.

5. Service / Maintenance

- » Use only soft clothes for cleaning. Do not use any detergents which may contain solvents or which scrub. Do not touch the lamp.
- » Keep the magnetic door strip clean.

Gloves

The gloves are wearing parts and have to be exchanged regularly. Due to the diversity of possible applications, a period cannot be defined. It is each user's obligation to check regularly the state of the gloves and to exchange them if necessary.

- » Remove the glove from the holding ring.

- » Slide the new glove onto the holding ring



Exchange of UV lamp

Generally, the lamp should be exchanged after approx. 10.000 operating hours. Use gloves in order not to destroy the lamp. The deinstallation / exchange should only be made by qualified persons.

Deinstallation of the UV lamp

- » Remove the UV installation from power supply / remove the power plug from the shockproof socket.
- » Loosen the fixing screw on the backside of the lamp housing by means of a fork wrench.
- » Loosen the cable connection and withdraw the cap.
- » Pull the UV lamp with its cable from both clip brackets.

- » Remove the lamp from the 4-PIN contact plug (Attention: may be hot!!) Only touch the UV radiator directly at the base (Risk of glass break!).

Installation of the UV lamp

- » Cleaning: Remove fingerprints on the lamp tube with an alcohol-infused soft cloth.
- » Put the UV radiator in the socket of the lamp connector cable of the ballast unit by means of the 4-PIN contact plug.
- » Clip the UV radiators in both lamp brackets.
- » Fix the lamp cable by means of the cable connection and realize strain relief
- » Tighten the safety screw for the cable by means of a fork wrench.
- » Put the installation into operation again. Reset the operating hours counter through pushing "Reset".

Note: Independent from the lamp age, the socket of the UV lamp seems brown under UV radiation effect. This does not affect the function.

6. Disturbances / Repair

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned glove box in its original packing to us. Then we will check whether a cost-effective repair is possible.

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7. Waste Disposal

The glove box as well as its packing are made of recyclable materials. The illuminant of the UV lamp contains mercury. Please waste the lamp at your local disposal point. Please take further information on the materials used from the following specifications.

8. Specifications

SICCO Glove Box Pure with Transfer Chamber and UV Sterilisation lamp

Cat. No.:	V 1740-08
Outerdimensions (W x H x D):	1210 x 700 x 600 mm
Weight:	ca. 20 kg
Usable interior space (W x H x D):	850 x 650 x 540
Usablevolume:	290 liters
Max. all-overload:	40 kg
Opening for gloves (O.D.):	Ø 190 mm
Distance between opening for gloves:	390 mm
Hinged lid opening (W x H):	850 x 425 mm
Opening for gas supply:	Ø 6 mm
Opening for gas outlet:	Ø 41 mm
Length of hose:	10 m tubing for gas-filling and 10 m aeration tubing

Details of transfer chamber

Outerdimensions (W x H x D):	320 x 330 x 260 mm
Usable interior space (W x H x D):	220 x 275 x 320 mm
Usablevolume:	24 Liter
Max.usableheight:	275 mm
Max. usablewidth:	200 mm
Distancetofloor:	40 mm
Max. all-overload:	3 kg

Details of UV lamp

The illuminant of the UV lamp contains mercury and must be wasted at local disposal points after exchange.

Electrical power:	36 W
Casing dimensions (L x W x H):	830 x 70 x 58 mm
Lamp power:	425 mA
Lampvoltage:	92 V
Lifetime:	10.000 h
Diameter:	15 mm
Length (withoutplug):	793 mm
Lengthplug:	8 mm
Irradiation width:	780 mm
Weight:	1,3 kg
Temperature max.:	40 °C
Heatingperiod:	30 seconds

9. Accessories

» see page 29-30

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CE Declaration of Conformity

EG-Declaration of Conformity in compliance with Machinery Directive 2006/42/EG

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herewith declares that the following product:

Glove Box Pure with Transfer Chamber, item number: V 1740-08

applies to the above mentioned directive including any amendments valid at the time this declaration was signed.

The low voltage directive 2006/95/EG has been fulfilled according to 2006/42/EG terms of protection objectives.



The following harmonized standards have been applied:

DIN EN ISO 12100	Safety of Machinery
ISO 12100:2010	Risk assessment and risk reduction
DIN EN 61326-1: 1997	Electrical Measuring, Control and Laboratory devices – EMV standards, Part 1: General requirements
DIN VDE 0100-100	Rules for the erection of high-voltage systems up to 1000 V
DIN EN61010-1:2010	Safety regulations for Electrical Measuring, Control and Laboratory devices – Part 1:General requirements
DIN EN 60335-1: 2012-10	Safety of electrical devices for domestic use and similar purposes – Part 1: General requirements

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Bohlender GmbH • Waltersberg 8 • 97947 Grünsfeld • Germany
BOHLENDER GmbH