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Operating Instructions

Mode d'emploi

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SICCO Gas Dosing Controller, Cat.-No.: V 1775-02

OPERATING INSTRUCTIONS

1. Supplied with	9
2. Security Advice	9
3. How it works	9
4. Taking into operation	10
5. Service / Maintenance	13
6. Disturbances / Repairs	13
7. Waste Disposal	14
8. Specifications	14
9. Accessories	14

Please read the information listed below thoroughly. We urgently ask you to observe the instructions about security, use and maintenance.

The SICCO Gas Dosing Controller controls the inlet of gas in dependence of the relative inside humidity of your glove box. The relative humidity is monitored by a sensor which has to be led into the glove box by using a cable lead-in.

1. Supplied with

- 1 Controller
- 1 Sensor cable
- 1 Power cable
- 1 Connection tubing 30 cm
- 1 Connection tubing 2,5 m
- 1 Operating instructions

2. Security Advice

- » Do not take a damaged Gas Dosing Controller into operation.
- » Keep the ventilation slots in the controller free.
- » The ambient temperature range of $-20\text{ }^{\circ}\text{C}$ to $+60\text{ }^{\circ}\text{C}$ shall not be exceeded.
- » Do not work with open door.
- » Just use the gas filling option in a well-vented room.
- » Transfer the extracted air in a controlled way, e.g. by using an extractor hood.
- » Make sure that the airflow of the extractor hood does not effect or ruin the process of gas filling.
- » Make sure to install all connecting kits for gas-filling.
- » Always use the smaller tubing as your tubing for gas-filling and the larger tubing as your aeration tubing.
- » The pressure for lead-in gas shall not exceed 0,3 bar.
- » Connect the aeration tubing at first, then the tubing for gas-filling. Disconnection has to be made in reverse order.
- » The cross section of the tubing must not be constricted.
- » Neither put the tubing for gas-filling nor the aeration tubing under tensile or compressive stress.
- » Avoid unanticipated reactions which might be caused by gas residues upon opening the swing-up door or the transfer chamber.
- » Always close the gas tap if regulation is no longer required.
- » These operating instructions cannot reflect all possible dangers which might occur during use. Users should check and observe all safety regulations.

3. How it works

The Gas Dosing Controller is installed between the gas tap and the glove box. Simply adjust the required relative humidity directly on the Gas Dosing Controller, the gas supply is open. When the adjusted value is reached, the supply on the Gas Dosing Controller closes. The sensor inside the work space constantly monitors humidity. Once humidity deviates for more than 1 % of the adjusted value, the gas supply automatically opens again after 3 minutes. Thus, humidity is kept on the same level. The time lag of 3 minutes prevents an unnecessary inlet of gas, for example, if the door of the glove box is only opened for a short period of time. For a deviation larger than 5 %, an alarm signal sounds. Such a deviation might occur, for example, if the door is permanently opened or if the gas bottle is empty.

Sicco Gas Dosing Controller

Front:

- » SET Button (10)
- » Alarm LED (11)
- » Temperature display (12)
- » Humidity display (13)
- » Up / Down Button (14)
- » On / Off Switch (15)

Back:

- » Gas Inlet (17)
- » Sensor cable port (18)
- » Gas Outlet (19)
- » Socket for power cable 230 V (20)


Front:

- » Flow meter

Back:


- » Inlet Flow meter
- » Outlet Flow meter

Info: All inlets and outlets as well as buttons not named are not necessary for the control of relative humidity.

4. Taking into operation

- » Connect the flow meter outlet (23) with the gas inlet (17). Push the 30 cm connection tubing (4) into the pneumatic quick coupling onto the bore ground.




- » **Tip:** To loosen tubing from pneumatic quick couplings, push the snap ring in direction of the housing. Pull the tubing out of the quick coupling while pushing the snap ring.



Sicco Gas Dosing Controller

- » For correct function, we recommend to use a SICCO Gas-saving valve, Cat.-No. V 1787-07.
- » Please note: gas which is lighter than air has to be lead in from above. Gas which is heavier than air from the bottom (density of air = 1,2041 kg/m).
- » Open the lower connector for gas-filling on the side panel of the work space through turning the sleeve nut (30) counter-clockwise. Remove the sealing disc (31) and keep the parts until they are reused.
- » To install the gas inlet tubing into the work space, uninstall a free cable lead-in on the back side of the glove box. Close the cable lead-in with the sleeve nut (30) and the sealing disc (31).
- » The cable lead-in is also available as an accessory, Cat.-No. V 1782-07 or V 1783-07 (Antistatic version).
- » Push the 2,5 m connection tubing for the gas inlet (5) with the parts removed above through the cable lead-in into the lower connector for gas-filling.
- » If required, cut the tubing with a sharp knife on a cut resistant pad.
- » Fix the four fastening screws by turning them carefully clockwise.
- » Push the tubing (5) into the gas outlet (19) onto the bore ground.


Sicco Gas Dosing Controller

- » Please note: In case of gas which is heavier than air (density of air = 1,2041 kg/m), make sure to lead in the sensor cable (2) into the upper cable lead-in so that the values are measured correctly.
- » If you use gas which is lighter than air, push the sensor cable (2) into the work space until the sensor cable measures from the lower part of the workspace.
- » Install the sensor cable (2) in the cable lead-in.


- » Fix the four fastening screws by turning them carefully clockwise.

- » Connect the sensor cable (2) with the sensor cable port (18). Please note that the plug has only one way to be plugged in due to its form.

- » Connect the flow meter inlet (22) with your gas tap by using an appropriate tubing (not indicated in the scope of delivery) with an O. D. of 8 mm. The gas tap is closed!
- » Adjust a pressure of between 0,1 and 0,3 mm on the pressure regulating valve.


- » Plug in the power cable (3) into the provided socket (20) to connect the device with a 230 V wall socket.

Sicco Gas Dosing Controller

- » Make sure that all openings are closed on your glove box and the gas is led out properly.
- » Open the gas tap.
- » Adjust the requested flow rate on the flow meter (21). The flow rate can be raised by turning counter clockwise. By turning clockwise, the flow rate is reduced.

- » **Info:** On switch-on of the Gas Dosing Controller, the sensor measures the current temperature and relative humidity inside the glove box. Those values are displayed. Depending on the adjusted relative humidity, the alarm signal possibly sounds immediately (if maybe the adjusted rH is very low). The alarm LED (11) flashes, an acoustic signal rings out. Reset the acoustic alarm through pushing the SET button (10). The alarm LED (11) continues to flash.
- » Switch on the Controller on the On / Off Switch (15).

- » Adjust the requested relative humidity. Push the UP-Button (14) as well as the Set-Button (10) at the same time. The currently adjusted relative humidity is shown on the left display.
- » Raise the value with the UP-Button(14), reduce the value with DOWN-Button (14).

- » Save the adjusted value through pushing again the Set-Button (10).
- » On the left display, the current temperature in the glove box is shown again, on the right display the current relative humidity. The controller starts to dose.

5. Service / Maintenance

Use only soft clothes for cleaning. Do not use any detergents which may contain solvents or which scrub. You can use all common neutral detergents (pH 7).

6. Disturbances / Repairs

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned Gas Dosing Controller in its original packing to us. Then we will check whether a cost-effective repair is possible.

Sicco Gas Dosing Controller

Objection	Reason and help
Relative humidity inside the work-space does not fall.	» The work space is not yet completely filled with gas. Info: If you use nitrogen in a SICCO Glove Box with Transfer Chamber, you can expect a reduction of rH of approx. 30% within 30 minutes (8 L/min. flow rate). » The flow meter (21) is closed or the adjusted value is too low for your application. Increase the flow rate. » The connection tubes are not connected properly. Check the tubing connections gas inlet (17), gas outlet (19), inlet flow meter (22) and outlet flow meter (21). » The adjusted value of rH is higher than the actually measured rH. Check the adjusted value by pushing the UP-Button (14) and Set-Button (10) at the same time. The current adjusted relative humidity is shown on the left display. Reduce the value with the DOWN Button (14). Save the value through. Save the adjusted value through pushing again the Set-Button (10). » Your Glove Box has leaks. Make sure that all openings are closed (ports, cable lead-ins). Close all openings. Check if all transfer chamber doors are closed. Check if the Swing-up door is closed tightly. If necessary, clean the door seals magnetic strips with a soft cloth.
The gas consumption is too high.	» Assemble the gas-saving valve V 1787-07 (accessories Cat.-No. V 1787-07) to the glove box and connect the outlet hose. The valve opens only to achieve a pressure compensation. This way the gas consumption can be reduced.

7. Waste Disposal

The Gas Dosing Controller as well as its packing are made of recyclable materials. Please take further information on the materials used from the following specifications:

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BOHLENDER GmbH