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Operating Instructions

Mode d'emploi

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SICCO Glove Boxes Air Stream, Cat.-No.: V 1745-08

OPERATING INSTRUCTIONS

1. Supplied with	17
2 . Security Advice	17
3. Taking into operation	18
3.1 Swing-up door	18
3.2 Cable lead-ins	19
3.3 Transfer Chamber	21
3.4 Iris-Ports	22
3.5 Gloves	22
4 Operation mode	23
4.1 Circulating air mode	23
4.2 Fresh air mode	24
4.3 Exhaust air mode	26
5. Handling of air filters	27
6. Service / Maintenance / Filter exchange	28
7. Disturbances / Repair	29
8. Waste Disposal	29
9. Specifications	29
10. Accessories	30

SICCO Glove Boxes Air Stream

Please read the information listed below thoroughly. We urgently ask you to observe the instructions about security, use and maintenance.

SICCO Glove Boxes provide safe handling of toxic materials, allow production of chemical compounds and protect sensitive substances from the environment. Optionally, the chamber can be flooded with an inert gas and a special atmosphere can be created. Models with side transfer chamber/s allow to introduce products without noticeable change of the atmosphere inside. All panels are made of high-quality and transparent plastic material: The aluminium profiles are joint together with connectors made of polyamide.

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- 1 SICCO Handschuhbox
- 1 Paar Handschuhe
- 1 HEPA H14 Filter (premounted)
- 3 Closing lids for supply and exhaust air (premounted)
- 1 Intake connector
- 1 Exhaust connector
- 1 Operating Instructions

2. Security Advice

- » Before changing into another operation mode, the glove box has to be disconnected from power supply.
- » Do not take a damaged glove box into operation.
- » Check before every use whether the glove box has been damaged.
- » Provide a safe stand for the glove box.
- » Do not work with opened swing-up door.
- » SICCO Glove Box shall be taken into operation only in rooms with sufficient ventilation.
- » Place the glove box in a reachable working height.
- » Do not expose the glove box to direct sunlight.
- » Do not store any solvents in the glove box.
- » Do not use any sharp or pointed tools inside the glove box.
- » Do not put any explosive or inflammable materials inside the glove box.
- » Do not put any radioactive products inside the glove box.
- » Transfer the extracted air in a controlled way, e.g. by using an extractor hood.
- » SICCO Glove Boxes are not suitable for overpressure or vacuum.
- » The ambient temperature range of $-20\text{ }^{\circ}\text{C}$ to $+70\text{ }^{\circ}\text{C}$ shall not be exceeded.
- » These operating instructions cannot reflect all possible dangers which might occur during use. Users should check and observe all safety regulations.

SICCO Glove Boxes Air Stream

3. Taking into operation

Remove the packing and the protection cover and place the glove box on a horizontal, plain bench.

Attention: Never lift the glove box at the transfer chamber.

Only lift the glove box at the bottom (see picture below):

LIFT HERE ONLY

LIFT HERE ONLY

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3.1 Swing-up door

» Always pull both handles to open the swing-up door. Swing the door panel up to the highest point until it stops.

Please note: Due to the strong magnetic lock it is harder to open the door the longer it was closed before. This is normal and can be compensated by pulling stronger than in everyday use.

» **Attention:** The door stop defines the maximum opening of the swing-up door. Put the door on the door stop. The adjoining picture shows the maximum possible opening of the door.

3.2 Cable lead-ins

The cable lead-ins on the back side of the chamber allow to lead in any standard cables including mounted plugs or any rigid tubing. Every cable or rigid tubing with an outer diameter from 0,1 mm up to 10 mm can be assembled regardless of the cross section the cable / tubing has.

» Unscrew the four cylinder screws. A suitable hexagon wrench is included.


» Remove the supporting ring (5) and the supporting disk with centring collar (6).

» Remove the silicone plug (7) and the sealing foil (8).

» Remove the clamp ring (9).



» Put your electrical device inside the glove box. Lead your cable through the supporting ring (5). Please pay attention to the supporting ring (5) being properly assembled with the notches in direction to yourself in order to fix the cylinder screws later.

SICCO Glove Boxes Air Stream

- » Guide the plug from inside the glove box through the housing (10) of the cable lead-in.
- » Pull the cable of your electrical device until the required position is reached. Slide the clamp ring (9) on the cable.
- » Afterwards slide the silicone plug (7) on the cable as well.
- » Then slide the supporting disk with centring collar (6) on the cable with the centring collar in direction to yourself.
- » Press the silicone plug (7) into the housing (10). While doing so put pressure to the housing (10) from outside as well.
- » Shove the supporting disk with centring collar (6) along the cable. Place the clamp ring with ridge on top of the silicone plug (7). Place the supporting ring (5) on the supporting disk with centring collar (6).

» Fix the four fastening screws by turning them carefully clockwise.

» The assembly of a rigid tubing is identical to the assembly of a cable. Just follow the above mentioned steps.

3.3 Transfer Chamber (only for models with transfer chamber).

The transfer chamber is used to introduce further products or tools during a running process. A contamination of the working space or a loss of the filled-in gas can be prevented.

» If the inner door of the transfer chamber is open, the interlock prevents accidental opening of the outer door.

» Put the products into the transfer chamber. A tray facilitates the removal of the items.
» Close the outer door.

» Open the inner door of the transfer chamber and remove the products or tools.

SICCO Glove Boxes Air Stream

- » Open the inner door of the transfer chamber and remove the products or tools.

3.4 Iris-Ports

- » The iris-ports can be replaced with new ones (available as accessory, cat.-no. V 1976-07). You can remove them by slightly rotating them and simultaneously pull them off the inlet port. Assemble the new iris-ports by slightly rotating them and simultaneously slide them on the inlet port.
- » Both inlet ports can be closed with a lid (available as accessory, cat.-no. V 1790-07) to protect the interior from small particles or dust. The lids do not close gas tight.
- » **Attention:** Do not close the inlet ports during exhaustion!

3.5 Gloves

The gloves are wearing parts and have to be exchanged regularly. Due to the diversity of possible applications, a period cannot be defined. It is each user's obligation to check regularly the state of the gloves and to exchange them if necessary.

- » Remove the glove from the holding ring.

- » Slide the new glove onto the holding ring.

4. Operation mode


In the Glove Box Air Stream, a HEPA H 14 Filter is assembled which restrains 99,995% of all particles larger than 0,1 μm . Depending on the operation mode, the circulating air, fresh air or exhaust air is filtered.

- » Wlf the rotary switch is pressed, the controlling ring on the switch shines and the circulation of air begins. Through turning the switch, the air flow quantity can be regulated.
- » We cannot recommend a reference point for the rotary controller as it depends from the application..
- » If you work with pulverized substances, the flow rate should be low to avoid turbulences.
- » You can see the resistance of the filter for air stream on the pressure gauge. If a value of approx. 600 Pa is reached, the filter is saturated and should be exchanged.

4.1 Circulating air mode

- » In this operation mode, the circulating air in the work space is continuously circulated. It is sucked on the back side, lead through the filter and blown back into the work space along the front panel (as delivered).

Before changing into another operation mode, the glove box has to be disconnected from power supply.

- » Before changing into another operation mode, the glove box has to be disconnected from power supply.

SICCO Glove Boxes Air Stream

» Open the lid (4) in the interior.


» Mount the gloves.

4.2 Fresh air mode

» In this operation mode, fresh air is sucked from the atmosphere, lead through the filter and blown back into the work space along the front panel. The interior air is displaced along the rear panel through incoming air and passed up to the atmosphere or to a connected exhaust system. Accessories for the connection to an existing exhaust connector are not included in the scope of delivery.

Before changing into another operation mode, the glove box has to be disconnected from power supply.

» For the fresh air operation mode, the closing lid (1) has to be removed. Take the closing lid on the grip and pull and turn it at the same time.

SICCO Glove Boxes Air Stream

» The angled intake connector has to be pushed into the bore through pulling and turning at the same time.

» The connector is correctly fixed if the red O-ring is no longer visible.

» The closing lid (2) must be removed, closing lid (3) is closed.

» Open the lid (4) in the interior.


» Mount the gloves.

4.3 Exhaust air mode

» In this operation mode, the motor sucks the interior air, leads it through the filter and passes the filtered air into the atmosphere through the exhaust opening. Ambient air is sucked through the hand openings so pressure compensation in the workspace is ensured. The flow movement to the interior avoids that air from the workspace can escape into the atmosphere.

Before changing into another operation mode, the glove box has to be disconnected from power supply.

» Suction opening (1) is closed with a lid.

» Remove the closing lids (2) and (3)

» The openings (2) and (3) are connected through the exhaust pipe.

» The pipe is correctly fixed if the red O-rings are no longer visible.

- » Close the lid (4) in the interior.

- » The gloves have to be removed for this operation mode. (You can obtain suitable Iris Port openings, Cat. No. V 1976-07 in the Accessories.)

5. Handling of air filters

There are some crucial principles to be considered when handling with air filters.

Storage

- » Filters have to be stored dry without having contact to walls and floors.
- » Filters should stay in their original packing until their use.
- » The storage places of air filters must not be able to provide a humid environment through condensing humidity.
- » Wet or humid filters should not be used at all and have to be disposed of. Even if the filters dried, they may be a health risk.

Handling of HEPA and compact filters

- » A Take care of not damaging the sealings during unpacking HEPA filters. Foamed sealings can be easily destroyed through shear forces.
- » Filters should stay in their original packing until their use.
- » Avoid touching the filtering surface. It is very sensitive can be damaged already through careless touching.
- » Even smallest damages of the filter surface make HEPA filters unusable.
- » HEPA filters are sometimes very heavy and bulky. Filter exchange and warehouse movements should probably be made by two people.
- » Check the tight fit of the filter in the filter body as canting can lead to bypasses.
- » For use of filters in Glove Boxes in safety-relevant sectors, we recommend to use a test device to guarantee a bypass-free installation (if a particle measuring is not necessary anyway).

Waste disposal of used filter components

- » Filter components from supply air operation have to be considered as household respectively commercial waste.
- » Exhaust and circulating air filters are normally polluted with hazardous filtrates and have to be disposed of as per these regulations.

6. Service / Maintenance / Filter exchange

- » Use only soft clothes for cleaning. Do not use any detergents which may contain solvents or which scrub.
- » Keep the magnetic door strip clean
- » We recommend to wear a protective overall, fine particle mask, protective glasses and gloves for a filter exchange.
- » The Glove Box has to be switched off during the filter exchange.
- » Make sure before the filter exchange that the filter components are free of processing errors.
- » Inspect the frames and sealings for wear.
- » Check if the installation shows traces of corrosion or other signs for condensing humidity.
- » If you find wet parts or corrosion, the problem has to be solved.
- » Used filter components are probably polluted with hazardous substances and/or concentrations of microbes and, thus, bear a risk.
- » Dispose of used filter components immediately. Do not store used filter components in areas where they can be mixed up with fresh filters or they can bear health risks (e. g. next to food).
- » Do not use used but still usable appearing filter components for other applications.

Before starting with work, read again all points on the handling of air filters attentively.

- » At first, interrupt power supply from the glove box, put it on the floor and open the inspection flap.

- » Sidewise to the filter, you will find two fasteners which have to be opened now.

- » Remove the filter through pulling it upwards out of the guiding.
- » Clean the surfaces of the filter body before installing a new filter, if necessary. Mount the new filter in reverse order.

7. Disturbances / Repair

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned glove box in its original packing to us. Then we will check whether a cost-effective repair is possible.

Beanstandung	Ursache und Abhilfe
The door does not close properly.	» There are particles like dust which are on the magnetic strip and / or on the door frame. Clean the magnetic strip and the frame with a soft cloth.
The gas consumption is too high.	» Assemble the gas-saving valve V 1787-07 (see accessories page 45) to the glove box and connect the outlet hose. The valve opens only to achieve a pressure compensation. This way the gas consumption can be reduced.
The glove box wobbles	» Place the glove box on a flat surface. Check that all rubber feet are present. If your work surface is uneven, the rubber feet can be shortened if necessary,so that the glove box is stable again.
The vacuum gauge shows more than 600 Pa.	» Turn the rotary controller lower. » Exchange the filter if necessary.

8. Waste Disposal

The glove box as well as its packing are made of recyclable materials. Please take further information on the materials used from the following specifications:

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Details of Transfer Chamber

Outer dimensions (W x H x D):	320 x 330 x 260 mm
Usable interior space (W x H x D):	220 x 275 x 320 mm
Usable volume:	24 litres
Max. usable height:	275 mm
Max. usable width:	200 mm
Distance to bottom:	40 mm
Max. all-over load:	3 kg

Details of the filter unit

Operation modes:	Circulating air, Fresh air, Exhaust air
Filter class as per EN 1822:209:	H14
Filtration performance (MPPS method):	99,995 % of all particles > 0,1 um are restrained
Dimensions of filter cartridge:	305 x 305 x 78 mm
Installation depth of filter cartridge:	78 - 83 mm

SICCO Glove Box with Transfer Chamber

Cat.-No.:	V 1745-08
Outer dimensions (W x H x D)	1210 x 700 x 600 mm
Weight:	ca. 42 kg
Usable interior space (W x H x D):	850 x 650 x 540
Usable volume:	290 litres
Max. all-over load:	40
Opening for gloves (O.D.):	Ø 190 mm
Distance between opening for gloves:	390 mm
Hinged lid opening (W x H):	850 x 425 mm
Material frame:	Aluminium with connectors made of polyamide
Material panels:	PMMA (polymethylmethacrylate)
Material base panel:	PE (polyethylene)

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Declaration of Conformity

EG-Declaration of Conformity in compliance with Machinery Directive 2006/42/EG

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herewith declares that the following product:

Glove Box Air Stream with Transfer Chamber, item number: V 1745-08

applies to the above mentioned directive including any amendments valid at the time this declaration was signed.

The low voltage directive 2006/95/EG has been fulfilled according to 2006/42/EG terms of protection objectives.

The following harmonized standards have been applied:

DIN EN ISO 12100	Safety of Machinery
ISO 12100:2010	Risk assessment and risk reduction
DIN VDE 0100-100	Rules for the erection of high-voltage systems up to 1000 V
DIN EN 61326-1: 2012	Electrical Measuring, Control and Laboratory devices – EMV standards
DIN EN 61010-1: 2010	Safety regulations for Electrical Measuring, Control and Laboratory devices – Part 1: General requirements
DIN EN 60335-1: 2012-10	Safety of electrical devices for domestic use and similar purposes – Part 1: General requirements

Restrictions:

Especially the following chapter of the a. m. harmonized standard DIN EN 60335-1: 2012-10 is partly not fulfilled by the installation component.

- **Protection against access of moving, active and isolating parts, even concerning air and creepage distance, (DIN EN 60335-1 Chap. 8, 20.2, 21.1, 21.2, 22.5)**

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As a manufacturing company we can make products to your specifications. Send us a drawing or sketch together with answers to the questions listed below.

*Return this form with contact details **by fax or e-mail** and we will be in touch to discuss details and give you a free, no obligation quote.*

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BOHLENDER GmbH