# Bedienungsanleitung Operating Instructions Mode d'emploi

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siehe Seite 37

|                                    |    |
|------------------------------------|----|
| <b>OPERATING INSTRUCTIONS</b>      |    |
| 1. SPECIFICATIONS.....             | 15 |
| 2. SCOPE OF DELIVERY .....         | 17 |
| 3. SAFETY INSTRUCTIONS.....        | 18 |
| 4. TAKING INTO OPERATION.....      | 18 |
| 5. MODE OF FUNCTION .....          | 20 |
| 6. SERVICE / MAINTENANCE .....     | 20 |
| 7. PROBLEMS / REPAIR .....         | 24 |
| 8. WASTE DISPOSAL.....             | 24 |
| 9. CUSTOMIZATIONS.....             | 25 |
| 10. SPARE PARTS / ACCESSORIES..... | 25 |

BOLA **Magnetic Stirrer Head**  
**P-MRK**  
Cat.-No.: C 520-..., 530-...

BOLA **Magnetic Stirrer Head**  
**P-MRK Flat Flange**  
Cat.-No.: C 522-..., 532

BOLA **Magnetic Stirrer Head**  
**P-MRK Rodaviss®**  
Cat.-No.: C 540-...

# BOLA Magnetic Stirrer Heads

BOLA Magnetic Stirrer Head can be used without additional lubricants and are suitable for long-term operation. Through the leading in a receiver, which is connected to the driving magnet, neither the guiding nor the stirrer shaft suffers from abrasion that could contaminate the medium. All spare parts can be disassembled for cleaning and be exchanged, if necessary.

Depending on the type of application, BOLA Magnetic Stirrer Heads can be obtained with ground joint or flat flange in two different types:

- »Type 1 is delivered with a receiver made of borosilicate glass. All parts coming into contact with the medium are metal-free as they are made of different chemically high-resistant plastic compounds or glass. This version is especially suitable for applications with chemically highly aggressive or highly pure products (Cat. No. C 520-..., C 522-...) .
- »Type 2 is delivered with a receiver made of Hastelloy®. Unlike glass, Hastelloy® is unbreakable. This allows to use stirrer shafts made of other materials such as stainless steel as well as a frequent disassembly of the stirrer shaft as there is no risk of glass breakage when leading in the shaft into the receiver (Cat. No. C 530-..., C 532-...) . BOLA Magnetic Stirrer Heads with Flat Flange can either be used with "QVF" flat flanges or with „SCHOTT" flat flanges. The medium is only exposed to metal-free materials which all have a very good chemical resistance.

## 1. Specifications » Bola Magnetic Stirrer Heads, Ground Joint Version

| Cat.-No. | For Ground Joint | Dimensions  | For shaft O.D | Material                                | Torque | Speed         | Connection to agitator |
|----------|------------------|-------------|---------------|-----------------------------------------|--------|---------------|------------------------|
| C 520-24 | NS 29/32         | Ø 50x243 mm | 8 mm          | PTFE, PTFE-compound, borosilicate glass | 90 Ncm | max. 1500 rpm | External square SW6    |
| C 520-28 |                  | Ø 50x241 mm | 10 mm         |                                         |        |               |                        |
| C 520-48 | NS 45/40         |             |               |                                         |        |               |                        |
| C 530-08 | NS 29/32         | Ø 50x243 mm | 8 mm          | PTFE, PTFE-compound, Hastelloy®         |        |               |                        |
| C 530-10 |                  | Ø 50x241 mm | 10 mm         |                                         |        |               |                        |
| C 530-15 | NS 45/40         |             |               |                                         |        |               |                        |

BOLA **Magnetic Stirrer Heads**

» **BOLA Magnetic Stirrer Heads, Flat Flange Version**

| Cat.-No. | For Ground Joint | Dimensions  | For shaft O.D | Material                                | Torque | Speed         | Connection to agitator |
|----------|------------------|-------------|---------------|-----------------------------------------|--------|---------------|------------------------|
| C 522-40 | DN 40            | Ø 62x241 mm | 10 mm         | PTFE, PTFE-compound, borosilicate glass | 90 Ncm | max. 1500 rpm | External square SW6    |
| C 522-50 | DN 50            | Ø 75x241 mm |               |                                         |        |               |                        |
| C 532-22 | DN 40            | Ø 62x241 mm | 10 mm         | PTFE, PTFE-compound, Hastelloy®         |        |               |                        |

» **BOLA Magnetic Stirrer Head (P-MRK) with Rodaviss®**

| Cat.-No. | For Ground Joint | Dimensions | For shaft O.D | Material                        | Torque | Speed       | Connection to agitator |
|----------|------------------|------------|---------------|---------------------------------|--------|-------------|------------------------|
| C 540-08 | NS 29/32         | Ø          | 8 mm          | PTFE, PTFE-compound, Hastelloy® | 90 Ncm | 1.500 1/min | External square SW6    |
| C 540-10 |                  | Ø          | 10 mm         |                                 |        |             |                        |

**2. Scope of Delivery**

- » 1 User Manuel
- » 1 BOLA Magnetic Stirrer Coupling (P-MRK)

**Component parts of Magnetic Stirrer Head for stirrer shaft dia. 8 mm.**

- » **Pos. 1** Driving Bell
- » **Pos. 2** Guide attachment
- » **Pos. 3** „SAFE-LAB“ or „Rodaviss“ nut  
(only for versions with ground joint)
- » **Pos. 4** Rotor with threaded pin
- » **Pos. 5** Thrust washer
- » **Pos. 6** Guide
- » **Pos. 7** Stirrer shaft receiver (made of borosilicate glass or Hastelloy®)
- » **Pos. 8** Clamping screw
- » **Pos. 9** O-Ring
- » **Pos. 10** Clamping wedge
- » **Pos. 11** Clamping nut

**Component parts of Magnetic Stirrer Head for stirrer shaft dia. 8 mm.**

Besides the above-mentioned items, the magnetic stirrer head for stirrer shaft dia. 8 mm includes the following component:

- » **Pos. 12** Reducing bush


See page 25 for a detailed spare parts list including catalogue numbers.

## BOLA **Magnetic Stirrer Heads**

### 3. Safety instructions

- » Do not take a damaged magnetic stirrer head into operation.
- » Only use the magnetic stirrer head up to the maximum speed of 1.500 rpm.
- » Never use the magnetic stirrer head (B) without globe stirrer coupling (A), Cat. No. C 399-12, or a comparable coupling.
- » There must be a gap of approx. 2 mm (see picture on the right side) between magnetic stirrer head (B) and globe stirrer coupling (A) for compensating heat expansion.
- » If the magnetic stirrer head falls down, there is a risk of injury during disassembly due to broken glass.
- » Only wear tight clothes to avoid being drawn in by rotating components.
- » If you have long hair, please wear a hairnet.
- » Please only use tight gloves.
- » Only use the magnetic stirrer head for products which are in a temperature range between - 200 °C to + 250 °C.

### 4. Taking into operation

- » Detach the BOLA „SAFE LAB“ Nut (3) by turning it counterclockwise (only ground joint version)

- » Center the „SAFE LAB“ nut (3) on the ground joint. The nut is providing additional safety, but it is not compulsory (only with Safe Lab).


- » We recommend to saw the upper part of BOLA Stirrer Shafts with a hacksaw. Remove the fin with a grinder or a file.


## BOLA Magnetic Stirrer Heads

- » Keep hold of the clamp screw (8) and open the clamp nut (11) by making counterclockwise turn.


- » Push the stirrer shaft through the clamp nut into the receiver. Attention: The receiver (7) consists of glass or Hastelloy®. Please be careful!

At least 40 mm of the PTFE-jacket must be pushed into the clamp nut in order to receive a safe guiding. But do not push the stirrer shaft as far as it will go, please leave a gap of approx. 5 mm.


- » Keep hold of the clamp screw (8) and tighten the clamp nut (11) clockwise. It is not necessary to use any tool for clamping the stirrer shaft.

### » **Brond joint version:**

Insert the assembled stirrer shaft through the ground joint into the lid of the reactor. The magnetic stirrer head should be placed on the nut. Keep hold of the magnetic stirrer head and tighten the nut clockwise until the magnetic stirrer head is placed free from backlash in the ground joint. If the stirrer shaft does not fit through the ground joint, push it through from below and fix it in the magnetic stirrer head.


## BOLA Magnetic Stirrer Heads

### Flat flange version:

- » Put a suitable gasket (20) (QVF or SCHOTT version) between the sealing surfaces. Put the magnetic stirrer head on the plane surface.


- » Assemble the flange connection according to manufacturer's instructions. In case of temperature fluctuation, we recommend to use a flat flange with adjusting springs (without illustration).

### Recommended for all magnetic stirrer heads

- » Place the globe stirrer coupling (A) between stirrer motor (C) and magnetic stirrer head (B).
- » Adjust the speed controller of your motor to 0 rpm. Approach to the requested speed by increasing it during operation. Please don't exceed the maximum speed of 1.500 rpm.

## 5. Mode of function

Due to the built-in permanent magnets in drive shaft and rotor, the torque is transmitted from the motor to the stirrer shaft in the reactor. The casing of the drive shaft does not turn..

## 6. Service / Maintenance


In case of strong contamination, the magnetic stirrer head should be disassembled and cleaned. PTFE and its compounds have a resistant, non-wetting surface which can be easily cleaned. Please do not use abrasive detergents since these might damage the surface and result in a milkiness of the surface. You can use all common neutral detergents (pH 7). For a stronger contamination we recommend to use an alkaline detergent up to pH 12. The BOLA Magnetic Stirrer Head can be used without additional lubricants.

## BOLA Magnetic Stirrer Heads


### How to disassemble the magnetic stirrer head

- » Keep hold of the clamp screw (8) and release the clamp nut (9) by making at least one counterclockwise rotation. You can also use two flat spanners size 22 for opening.

**Hint:** A jammed stirrer shaft can be released by swaying it carefully.


- » If the stirrer shaft can't be released, the clamp nut (9) should be completely screwed off. Attention: O-ring and clamp ring can easily get lost.


- » Keep hold of the drive shaft (1) by holding the casing and the square of the drive shaft. Completely turn off the guiding (6) clockwise. Attention: left-hand thread.

- » After nine rotations, the complete rotor assembly can be pulled out.


- » Screw off the setscrew in the rotor counterclockwise. Use a screwdriver with a width of 4.5 mm.


## BOLA **Magnetic Stirrer Heads**


- » Detach receiver (7) with guiding (6) from the rotor (4).


- » Detach the wearing disc (5) from the rotor (4) by using your fingernail.


- » Remove the guiding (6) from the receiver (7).

- » Pull the clamp screw (8) from the receiver (7). Here you will feel a slight resistance.


- » Remove the guiding assembly (2) from the drive shaft (1) by turning it counterclockwise. You can also use a flat spanner size 41 (Cat. No. C 520-48: flat spanner size 56).

## BOLA **Magnetic Stirrer Heads**

**After cleaning, assembly must be made in reverse order. Particular attention should be paid to following issues:**


- » Push the clamp screw (8) on the receiver (7) until it clicks into place.


- » Put the wearing ring (5) on the rotor (4). The hexagonal bolt on the wearing ring must fit into the hexagon of the rotor. You do not need any tools for assembly.


- » Center the ground surface of the receiver (7) into the threaded bore of the rotor (4).

- » Screw the setscrew with a screwdriver (width 4.5 mm) into the rotor until it rests on the surface. Turn the screwdriver clockwise. The setscrew must now be completely sunk in the rotor. Attention: Do not damage the setscrew by turning it too tight!

- » The guiding (6) must now be easily turnable and must not be clamped between wearing ring (5) and clamp screw (8).

**7. Problems / Repair**

In case of problems you cannot handle, please do not hesitate to contact us by phone. Of course we will try to help you and give advice to find a solution. If the problem might not be solved this way, the cleaned and originally packed magnetic stirrer head can be shipped to us. We will check then if a cost-effective repair is possible.

| Problem                                     | Reason and repair                                                                                                                                                                                                                                                                  |
|---------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The magnetic stirrer head is rough-running. | Dirt particles might have entered the magnetic stirrer head. Detach and clean the magnetic stirrer head.                                                                                                                                                                           |
|                                             | The receiver (7) made of borosilicate glass might be broken. Detach the magnetic stirrer head and replace the receiver.<br><b>Attention: Risk of injury because of broken glass!</b> Check all components for subsequent damage from broken glass. Replace all damaged components. |

**8. Waste disposal**

The packing as well as the magnetic stirrer head are made of recyclable materials. Please take further details about the materials from the following specifications.

### **9. Customizations**

- » If you need a different version of the BOLA Magnetic Stirrer Head, do not hesitate to send us your inquiry including all requested data. As a manufacturer, we are able to offer an individual custom manufacture. In general, following modifications of the standard magnetic stirrer head are possible:
- » A smaller receiver.
- » A different type of drive.
- » A different flat flange size.

Our standard stirrer shafts can be modified in almost every feature (stirring dia., length, shaft dia., shape of blade ....).

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herewith declares that the following products:

**Device type:** Magnetic Stirrer Couplings and Magnetic Stirrer Heads

**Model designation:**

Magnetic Stirrer Couplings (P-MRK), Article Group: C 520- ...

Magnetic Stirrer Couplings (P-MRK) Flat Flange, Article Group: C 522- ...

Magnetic Stirrer Heads with Joint (MRK), Article Group: C 450- ...

Magnetic Stirrer Heads with Flange (MRK), Article Group: C 454- ...

apply to the above mentioned directive including any amendments valid at the time this declaration was signed.

The following directives have been applied:

**Directive 2014/30/EU**

**Directive relating to electromagnetic compatibility, General requirements**

The low voltage directive 2006/95/EC has been fulfilled according to 2006/42/EC terms of protection objectives.

The following harmonized standards have been applied:

DIN EN ISO 12100

Safety of machinery

ISO 12100:2010

Risk assessment and risk reduction

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Hastelloy® is a registered trade mark of Haynes International.

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Rodaviss® is a registered trade mark of V.E.R.A.L.

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BOHLENDER GmbH