Bedienungsanleitung
Operating Instructions
Mode d'emploi

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BOLA Distributors for Reaction Vessels

BOLA Distributors for Reaction Vessels (R), Cat.-No.: D 744-16 and D 744-24.

OPERATING INSTRUCTIONS

1. Supplied with	8
2. Security instructions	8
3. Taking into operation	8
4. Mode of Function	9
5. Service / Maintenance	10
6. Disturbances / Repairs	11
7. Waste disposal	11
8. Specifications	11
9. Spare Parts / Accessories	12
10. Custom manufacture	12

Please read the information listed below thoroughly. We urgently ask you to observe the information about security, use and maintenance.


BOLA Distributors for Reaction Vessels (R) allow to use a standard laboratory flask as a reaction vessel. They are suitable for all DURAN glass bottles with thread GLS 80.

BOLA Distributors for Reaction Vessels (R) are available with an integrated stirrer bearing for either shaft dia. 8 mm or 10 mm.

BOLA Distributors for Reaction Vessels

1. Supplied with

- 1 Distributor for Reaction Vessels (1)
- 1 Shaft Guide (2)
- 1 Gasket (3)
- 1 Screw Cap (4)

2. Security instructions

This operating instruction cannot give advice for all problems that may occur. Users should check the suitability and compliance of parts and materials with security and health instructions before use.


- » Do not take a damaged Distributor for Reaction Vessels (R) into operation.
- » Any product inside the vessel shall not exceed the operating temperatures of -20 °C up to + 140 °C.
- » Check whether any inserted probes or tubes may collide with the stirrer shaft before starting the agitator.
- » Start your agitator at low speed and accelerate slowly up to the requested speed.
- » Use only suitable laboratory screw joints (e. g. BOLA Cat.-No. D 629-.. und D 630-..) to connect probes or tubing to the GL sockets.
- » Close unused GL sockets with a suitable GL screw cap (e. g. BOLA Cat.-No. H 993-18 or H 994-18).
- » Check whether the lid is fixed tightly on your vessel before starting the agitator.
- » Provide a safe stand for your vessel.
- » After exchange of the shaft guide make sure that you fix the shaft with the suitable screw cap.
- » Use for stirrers with dia. 8 mm the BOLA Screw Cap C 425-84 and for stirrers with dia. 10 mm the screw cap C 425-86.

3. Taking into operation

- » Loosen the screw cap GL 25 on the central socket by turning it one turn to the left.

- » Push the shaft through the shaft guide. Once you have reached the requested height, tighten the GL 25 Screw Cap. The shaft is now fixed and sealed.

BOLA Distributors for Reaction Vessels

- » Further components can be connected to the lateral GL sockets by using suitable BOLA Laboratory Screw Joints for thread GL 18, e.g. cat.-no. D629-xx. (Screw joints are not included in delivery). BOLA Laboratory Screw Joints consist of a screw cap (5), a clamp ring (6), a tapered ring (7) and a sealing ring (8).

4. Mode of Function

» Lateral Sockets GL 18

The lateral threaded sockets can be used for insertion of tubes, tubing or probes with a max. OD of 10 mm. Please use suitable laboratory screw joints e.g. BOLA Cat.-No. D 629-xx (not included in delivery).

» Central Socket GL 25

A shaft guide can be put into the central socket and thus a stirrer shaft can be inserted, fixed and sealed. For shaft dia. 8 mm please use the shaft guide C 425-58 and screw cap C 425-84.

For shaft dia. 10 mm please use the shaft guide C 425-59 and screw cap C 425-86.

- » The central socket GL 25 is used for insertion of a stirrer shaft.

- » A BOLA Multiple Distributor cat.-no. D 619-08 can be fixed on the central socket. Tubes and tubing can be connect to in total 4 sockets GL 18 and 3 sockets GL 14 (distributor not included in delivery).

- » By means of a BOLA Laboratory Screw Joint GL 25 (cat.-no. D 630-xx), hard-walled tubes and tubing or probes etc. can be connected to the central neck (laboratory screw joint not included in delivery).

This is only a small choice of possible connections and enlargements with the modular BOL GL Screw Joint System.

BOLA Distributors for Reaction Vessels

5. Service / Maintenance

Do not use any abrasive detergents as these can damage the surface. All common neutral cleaning agents (pH 7) can be used. For heavy contaminations, we recommend to use an alkaline detergent (up to pH 12).

Due to friction of the stirrer shaft during permanent use, the shaft guide as well as the gasket become worn and lose their function. With regards to the high number of possible applications (speed, viscosity and corrosiveness of the product, temperature etc.), we are not able to advice a maintenance rate. Please check regularly the state of the components. It is the user's responsibility to check and exchange them regularly.

The shaft guide and the gasket are wearing parts. They can be exchanged as follows:


- » Loosen the screw cap GL 25 on the central socket by turning it to the left.


- » Push the shaft guide (2) from the lower side of the lid (1).

- » Then pull the shaft guide (2) out of the lid (1).


- » Remove the gasket (3) from the shaft guide (2). Please use tweezers to pull the gasket out.

BOLA Distributors for Reaction Vessels

» Put a new gasket (3) into the shaft guide (2). The black O-ring has to be visible on the top.

» Push the shaft guide (2) back into the lid (1).

» Fix the screw cap on the central socket by turning it to the right.

6. Disturbances / Repairs

If you are not able to handle and to eliminate a disturbance by yourself, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned unit in the original packing to us. Then we will check whether a cost-effective repair is possible.

7. Waste disposal

BOLA Distributors for Reaction Vessels® as well as their packing are made of recyclable materials. Please take further information on the materials used from the following specifications:

8. Specifications

BOLA Distributors for Reaction Vessels (R)

Cat.-No.	Guide for shaft O.D.	For bottle thread	Materials	Temperature Resistance	Sockets
D 744-16	8 mm	GLS 80	PTFE, PP, PPS, FPM (Viton)	- 20 °C up to + 140 °C	1 x GL 25
D 744-24	10 mm				4 x GL 18

BOLA Distributors for Reaction Vessels

9. Spare Parts / Accessories

» see page 18

10. Custom manufacture

The BOLA Distributors for Reaction Vessels (R) can also be offered in custom sizes and versions. Simply send us your enquiry with the requested design. As manufacturer we are able to answer enquiries for custom-made items.

In general the following adaptations of BOLA Distributors for Reaction Vessels (R) are possible:

- » Different bottle thread.
- » Different numbers and sizes of central or lateral GL necks.
- » Different material.

Please contact us if you have questions regarding a custom manufacture. We will be pleased to help.

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BOHLENDER GmbH