Bedienungsanleitung Operating Instructions Mode d'emploi

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OPERATING INSTRUCTIONS

1. SUPPLIED WITH	8
2. SECURITY INSTRUCTIONS	8
3. TAKING INTO OPERATION	8
4. MODE OF FUNCTION	8
5. SERVICE / MAINTENANCE	9
6. DISTURBANCES / REPAIRS.....	9
7. WASTE DISPOSAL	9
8. SPECIFICATIONS.....	10
9. CUSTOMIZED TEMPERATURE PROBES.....	11

BOLA **Temperature Probes PT 1000**
Cat-No.: P 1950-...

BOLA **Immersion Probes PT 1000**
Cat-No.: P 1980-20

BOLA **Temperature Probes 1000**
LEMO® Cat-No.: P 1962-...

BOLA **Immersion Probes PT 1000**
LEMO® Cat-No.: P 1992-20

Please read the information listed below thoroughly. We urgently ask you to observe the information about security and use.

BOLA Temperature Probes are PT 1000 probes. They are ideal for measuring temperatures in aggressive liquids. They are available with connectors type LEMO®, plug size 1. The LEMO® connectors are already mounted and can be connected to every suitable measuring instrument.

BOLA Temperature Probes PT 100

1. Supplied with

- 1 Temperature Probe PT 1000 with cable length 1.500 mm or 1 Immersion Probe with cable length 4.000mm (both alternatively with or without LEMO® plug)
- 1 Operating Instructions


2. Security instructions

The temperature probes are designed to make measurements in aggressive liquids. The operating instructions cannot show all security problems which can arise eventually. It is up to the user to make sure that all safety and health regulations are kept and that the corresponding restrictions are stipulated before use.

- » Please do not take a damaged temperature probe into operation.
- » Please take the response times of the temperature probe into consideration during use. Typical response times T 50 and T 90 (50% respectively 90% of the final measuring value) are to be found under point 8 Specifications on page 10
- » Please do not exceed the operating temperature range of maximum +250°C and minimum -200°C.
- » Please verify that the LEMO® connectors are linked correctly to the measuring instrument.

3. Taking into operation

- » The temperature / immersion probes only have to be connected to a suitable measuring instrument. After that, you can start measuring.
- » For temperature and immersion probes without LEMO® connectors, the following connecting plan has to be considered. If you want to mount your own connector, please consider the correct position of the pins.

4. Mode of function

BOLA temperature probes are sensors based on platinum resistance changes under temperature influence. At a temperature of 0°C they have an electrical resistance of 1000 ohm. Therefore the additional resistance caused by the measuring line plays a subsidiary role. The temperature probes PT 1000 are manufactured in conformance to class A. The deviation of the measured temperature is equal to or less than $\pm 0,15^{\circ}\text{C} + (0,002 \times T)$. Typical deviations are $\pm 0,15^{\circ}\text{C}$ at a temperature of 0° C and $\pm 0,35^{\circ}\text{C}$ at a temperature of 100° C.

The measuring sensor itself is located at the end of a PTFE-encapsulated stainless steel tube (material code 1.4571). This tube provides certain rigidity, but can manually be bent to the requested form. The seamless PTFE-encapsulation provides an almost universal chemical resistance.

» BOLA Temperature probe PT 1000

The probes have tapered tips which reduce the response time considerably. The cable is encapsulated with PFA and connected tightly to the temperature probe. The 2-wire cable has a length of 1.500 mm and is equipped with a LEMO® connector, plug size 1.

» BOLA Immersion probe PT 1000

This probe can be immersed completely in aggressive liquids and gases. The 2-wire cable is coated with white PFA. The cable and the probe are resistant to virtually all chemicals. In addition, the total immersion probe can be used for measurement of temperature in autoclaves. They are equipped with a LEMO® connector, plug size 1 and 2-wire cable with a length of 4.000 mm.

5. Service / Maintenance

In case of a strong contamination, the temperature probe should be cleaned after use. Generally, all fluoroplastic materials, PTFE, PFA and FEP, have smooth, non-wetting surfaces and can usually be cleaned without any problems. Please do not use any abrasive scouring agents as they might damage the surfaces. You may use all neutral detergents (pH 7). For a stronger contamination we recommend the use of an alkaline detergent up to pH 12. Apart from that, the probes are maintenance-free.

6. Disturbances / Repairs

If the probe should not work properly, please check first if the LEMO® connector is linked correctly to the measuring instrument. Also a damage of the connecting cable or of the sensor can lead to disturbances. Should there be a disturbance you cannot handle, please do not hesitate to call us. Of course we will try to help you this way. If we do not find a solution to your problem then, we will ask you to return the cleaned probe in the original packing to us. We will check if a cost-effective repair is possible.

7. Waste disposal

The probes as well as their packing are made of recyclable materials. Please take further information about the materials used from the following specifications:

BOLA Temperature Probes PT 100

8. Specifications

» BOLA Temperature Probes PT 1000

Cat.-No.:	P 1950-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA- coated cable
Temperature range	max. +250°C, min. -50°C
Resistance	1000 Ohm at 0°C
Probe-Ø	8 mm
Useful length	200 / 300 / 500 / 600 mm
Total length	approx. 260 / 360 / 560 / 660 mm
Tip-Ø	5,5 mm
Collar ring-Ø	approx. 12 mm
Cable length	1.500 mm
Response time	T 50: 16 - 18 s ; T 90: 47 - 50 s

» BOLA Temperature Probes PT 1000 LEMO® with plug

BOLA Cat.-No.:	P 1962-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA- coated cable Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	1000 Ohm at 0°C
Probe-Ø	8 mm
Useful length	200 / 300 / 500 / 600 mm
Total length	approx. 260 / 360 / 560 / 660 mm
Tip-Ø	5,5 mm
Collar ring-Ø	approx. 12 mm
Cable length	1.500 mm
Response time	T 50: 16 - 18 s ; T 90: 47 - 50 s

» BOLA Temperature Probes PT 1000 LEMO® with plug

BOLA Cat.-No.:	P 1980-20
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA- coated cable
Temperature range	max. +250°C, min. -50°C
Resistance	1000 Ohm at 0°C
Probe-Ø	4 mm
Probe length	approx. 50 mm
Cable length	4.000 mm
Response time	T 50: 16 - 18 s ; T 90: 47 - 50 s

BOLA Temperature Probes PT 100

» BOLA Immersion Probes PT 1000 LEMO® with plug	
BOLA Cat.-No.:	P 1992-20
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA- coated cable Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	1000 Ohm at 0°C
Probe-Ø	4 mm
Probe length	approx. 50 mm
Cable length	4.000 mm
Response time	T 50: 16 - 18 s ; T 90: 47 - 50 s

9. Customized Temperature Probes

Should you be in need of a temperature or immersion probe with other dimensions or connectors, simply send us an inquiry with the requested specification. As manufacturers we are the best contact for custom production.

Generally, the following modifications of our standard temperature probes can be offered:

- » Special dimensions of the thermo element (different probe diameter, useful or total length).
- » Different cable lengths..
- » Different kind of LEMO® connectors – a connector of your choice can be mounted.
- » We can coat your temperature probes with a PTFE heat shrinkable tubing so that they have the chemical resistance of PTFE. Because of the thin coating, the probe or thermometer has slower response behaviour. For coating, your probe has to withstand a short-time temperature of minimum +250 °C.

Please do not hesitate to contact us if you have any questions regarding custom manufacture.

» Registered Trademarks

LEMO® is a registered trademark by Lemo.

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Certificate

CE Declaration of Conformity

in compliance with Machinery Directive 2006/42/EG

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Herewith declares that the following products:

Type designation: Temperature Probes Lemo Compact, product group: P 1730- ...
Double Temperature Probes Lemo Compact, product group: P 1740- ...
Temperature Probes, product group: P 1750- ...
Temperature Probes Lemo, product group: P 1760- ... / P 1762- ...
Total Immersion Probes, product group: P 1780- ...
Total Immersion Probes Lemo, product group: P 1790- ... / P 1792- ...
Temperature Probes PT 1000, product group: P 1950- ...
Temperature Probes PT 1000 Lemo, product group: P 1962- ...
Total Immersion Probes PT 1000, product group: P 1980- ...
Total Immersion Probes PT 1000 Lemo, product group: P 1992- ...

apply to the above mentioned directive including any amendments valid at the time this declaration was signed.

The following directives have been applied:

Directive 2014/30 EU

**Directive relating to electromagnetic compatibility, General requirements,
The low voltage directive 2006/95/EC has been fulfilled according to 2006/42/EC terms
of protection objectives.**

The following harmonized standards have been applied:

DIN EN ISO 12100 – Safety of machinery

ISO 12100:2010 – Risk assessment and risk reduction

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BOHLENDER GmbH