Bedienungsanleitung
Operating Instructions
Mode d'emploi

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» siehe Seite 18

SICCO Vacuum-Desiccators, Cat.-No.: V 1880-04, V 1880-08 and V 1880-12

OPERATING INSTRUCTIONS

1. Content of delivery	8
2. Security instructions	8
3. Taking into operation	8
4. Mode of function	10
5. Service / Maintenance	10
6. Disturbances / Repairs	10
7. Waste disposal	11
8. Specifications	12
9. Spare/special parts	12

SICCO Vacuum-Desiccator
Cat.-No.: V 1880-04

SICCO Vacuum-Desiccator
Cat.-No.: V 1880-08

SICCO Vacuum-Desiccator
Cat.-No.: V 1880-12

Please read the information listed below thoroughly. We urgently ask you to observe the information about security, use and maintenance.

SICCO Vacuum-Desiccators are suitable for storing sensitive products in a vacuum. The vacuum inside the desiccator ensures rapid drying of the stored products. These systems allow storage in an oxygen and carbon dioxide free environment and are, for example, suitable for storage of semiconductor samples.

1. Content of delivery

- 1 SICCO Vacuum Desiccator
- shelves (V 1880-04: 2 shelves / V 1880-08: 3 shelves / V 1880-12: 4 shelves)
- 1 vacuum gauge, mounted
- 2 needle valves with hose connector, mounted
- 1 Operating instruction

2. Security instructions

- » Please do not take a damaged desiccator into operation.
- » Provide a safe stand for the desiccator.
- » Place the desiccator in a reachable working height.
- » Do not put any explosive or inflammable products into the desiccator.
- » Do not store any solvents or acids in the desiccator.
- » SICCO Desiccators are not suitable for overpressure
- » Only evacuate the desiccator at room temperature.
- » Don't produce more than 10-4 Torr negative pressure.
- » Never open the door when the desiccator is evacuated.
- » At normal pressure, the desiccator can be used in ambient temperatures between -20°C and $+70^{\circ}\text{C}$.
- » Do not expose the desiccator to direct sunlight.
- » Remove only one shelf at once.

3. Taking into operation

- » Remove the protection cover.



- » Tighten the hose connectors (1) and (2) carefully by turning clockwise with a flat wrench size 15.
- » **Attention:** You must have a sufficient counterholding force to prevent damages at the bonding surface of the valve.
 - 1: Valve for evacuation
 - 2: Valve for ventilation
 - 3: Vacuum gauge

- » Push a hose clamp (not included in delivery) for approx. 10 cm on your hose.

- » Push the hose over the connector of the valve (1).


Attention: Never turn the hose in counterclockwise direction while pushing it over the connector.

- » **Tip:** The hose can be pushed easier if it is heated before (for example with a hair dryer or with warm water).


Attention: Never turn the hose in counterclockwise direction for removing.

- » Fix the hose with the hose clamp.
- » Connect the hose to the vacuum pump according to the operating instructions of the pump.


We recommend to integrate a cold trap between desiccator and pump to protect the pump from condensates and aggressive products


- » Open the valve (1) which is connected to the vacuum pump by turning completely counterclockwise.
- » Close the other valve (2) by turning carefully clockwise until you feel a slight resistance.

- » Put your components into the vacuum desiccator.
- » Close the door and block all latches (4) by turning counterclockwise.

- » Activate your vacuum pump and monitor the evacuation procedure. The mounted vacuum gauge (3) indicates the inside pressure.
- » Close the valve (1) by turning clockwise and switch off your vacuum pump as soon as the requested pressure is reached.


- » **Attention:** The position of the latches (4) must not be changed as long as the desiccator is under vacuum. Otherwise they deform during pressure compensation and lose function.

After drying, the stored components can be taken out as follows:


- » Open the valve for ventilation (2) carefully. The valve is cooled during ventilation. This does not affect function.

- » Open all latches (4) by turning clockwise.

Attention: the door can only be opened after pressure compensation is completed.

- » Open the door and take out the stored products.

4. Mode of function

- » Due to the connection to a vacuum pump, the inside pressure of the desiccator is decreased and drying is speeded up.
- » The vacuum gauge indicates the inside pressure.
- » The door cannot be opened under vacuum.
- » To prevent the air humidity getting into the desiccator when putting in or taking out parts, the door should be opened only for a very short period of time.

5. Service / Maintenance


- » Only use soft clothes and no detergents which contain solvents or which scrub.
- » Only use the desiccator in clean and dustfree conditions. Never expose it to direct insulation or other sources of heat.

6. Disturbances / Repairs

The following information to clear possible disturbances shall help you to recognize or to eliminate possible reasons for the disturbances by yourselves. Should there be a disturbance you cannot handle, please do not hesitate to contact us. Of course we will try to help you in this way. If then we do not find a solution to your problem, we ask you to return the cleaned desiccator into the original packing to us. Then we will check if a cost-effective repair is possible.

Objection	Reason and help
The air pressure inside the desiccator does not decrease.	» The door is not closed tightly. Close the door and lock all latches. » Check the circular silicone seal for damages. If necessary, renew the silicone seal. » The valve for evacuation (1) which is connected to the pump is not open. Open the valve (1). » The valve for ventilation (2) is still open. Close the valve (2). » Check your connectors, hose and pump for leakages. Seal the leakages. Renew the hose.
The vacuum gauge is defective	» Exchange it as follows: You can use all vacuum gauges with G 1/4" threads)

- » Remove the defective vacuum gauge by turning counterclockwise using a fork wrench size 14.
Attention: Do not damage the vacuum gauge by turning too strong.

- » Apply a suitable adhesive on the thread.
 - » Screw the new vacuum gauge into the desiccator by turning clockwise using a fork wrench size 14.
- Attention:** Do not damage the vacuum gauge by turning too strong.

7. Waste disposal

The SICCO vacuum desiccators as well as the packing consist of recyclable materials. Please take further information about the materials used from the following specifications.

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BOHLENDER GmbH