

b.safe SYSTEM



New in our range:
b.safe Funnel Universal
R 803-04

b.safe SYSTEM

Hinged Lid

- Protection from harmful vapors

Funnel

- made of PE-EX, for disposal of liquid waste

Collective Connection

- for capillaries and tubing of up to 4 HPLCS

Screw Cap GL45

- for connection to a **b.safe** wall or desk-top mount



We manufacture waste disposal solutions for all of your challenges.

Big advantages put in small details:



Ask for your free b.safe SYSTEM catalogue

Scan QR-Code and
order your catalogue
www.bsafede.de/catalogue-request

