

b.safe SYSTEM



b.safe Mini-Display Leak

Cat. No.: R 658-21, R 658-22, R 658-24

b.safe Leak-Sensor

Cat. No.: R 637-01

NEW
in our
assortment

b.safe SYSTEM

Fill level monitoring 2.0: With innovative leakage protection


- Safety without compromises
- Leakage? No chance!
- No undetected leakage of solvents



- The **Leak-Sensor** triggers the alarm on the Mini-Display in the event of a leakage of a collecting container
- **Full canister:** sensor LED flashes red, a pulsating alarm tone sounds
- **Leaked liquid in the collecting tray:** sensor LED lights up red, a continuous alarm tone sounds



Have everything under control!
Our Leak-Sensor is the perfect supplement for collecting trays in safety cabinets.



We manufacture disposal solutions for your challenges



b.safe Mini-Display Leak – Monitors not only the level in connected collection containers but also the collecting tray in the safety cabinet. For this purpose, an additional **b.safe Leak-Sensor** is connected

b.safe Mini-Display Leak, 1-channel

» Cat. No. R 658-21

b.safe Mini-Display Leak, 2-channel

» Cat. No. R 658-22

b.safe Mini-Display Leak Splitter 2-channel

» Cat. No. R 658-24



The leak sensor triggers the alarm on the Mini-display in the event of a leakage of a collecting container

b.safe Leak-Sensor

» Cat. No. R 637-01

Ask for your free **b.safe SYSTEM** catalogue

Scan QR-Code and
order your catalogue

www.bsafe.de/catalogue-request

